



WONDFO BIOTECH
WeAreWorkingForYour **Health**

Distributed by:  MEDIVEN
MEDICAL INNOVATION VENTURES

2019-nCoV Antigen Test (Lateral Flow Method)

HOME NASAL SWAB SELF-TESTING



Guangzhou Wondfo Biotech Co., Ltd.
Add: No. 8 Lizhishan Road, Science City, Luogang District,
510663, Guangzhou, P.R. China
Tel: +86-20-32296083 Fax: +86-20-32296063
E-mail: sales@wondfo.com.cn Websites: en.wondfo.com.cn
Wondfo reserves all rights. All contents of this brochure are copyrighted by Wondfo.



WONDFO SELF-TESTING

Your best choice to monitor infection status for COVID-19.

- CE certified for home antigen self-test!
- Simple nasal swab sampling.
- Easy to use, give instant and reliable result after 15 mins.
- Various specifications available to meet different needs.

THE IMPORTANCE OF SELF-TESTING



Home self-testing is emerging as routine screening measure for pandemic control. Most of the EU countries have authorized the home antigen test.





Real-time monitoring of infection status for self and family members can better protect the community.




Self-performed test reduce the burden of healthcare facility (e.g. diagnostic center, central laboratory).


TEST KIT COMPONENTS



TEST PROCEDURES



Result interpretation



PERFORMANCE

Reagent	Confirmed by Referenced method PCR		Total
	Positive	Negative	
2019-nCoV Antigen Test (Lateral Flow Method)	Positive	230	231
	Negative	21	643
Total	251	623	874

CT Value	CT ≤25	CT ≤30	CT ≤33
Sensitivity	100%	98.92%	92.61%

Test's Sensitivity: 91.63%
Test's Specificity: 99.84%
Total agreement: 97.48%

ORDER INFORMATION

Catalog No.	Packing Size	Reaction Time	Shelf Life	Sample Type	Storage Temperature	Qualification
W634P0025~W634P0028	2/5/10/20 T/kit	15~20 minutes	24 months	Nasal swab	2~30°C	CE certified